Infusion/Withdrawal Pumps: KDS 310 Plus



Tomopal Part #	Manufacturer #	Voltage
781-3113 USA	781311 USA	115V~, 0.1A
789-3113 CE	789311 CE	230V~,0.05A

Features:

- Nanoliter Syringe Pump 1 μl-250 μl
- Designed to hold glass syringes from 1 to 100 microliters.
- Dispense Volume Mode: Pump keeps track of the volume dispensed and automatically stops the pump when a set target volume is reached.
- Run Mode: Pump runs at the set flow rate until manually stopped.
- Footswitch interface for hands free operation (optional).
- Min. Flow Rate (10 μ l syringe): 0.001 μ l/min
- Max. Flow Rate (10 µl syringe): 145.6 ml/min

Tel: 916.429.7240 Fax: 916.429.2701