

Infusion Pumps: KDS 100



| Tomopal Part # | Manufacturer # | Voltage |
|----------------|----------------|-------------|
| 781-1003 USA | 781100 USA | 115V~, 0.1A |
| 789-1003 CE | 789100 CE | 230V~,0.05A |

Features:

- Single-Syringe Infusion Pump
- \bullet The KDS 100 is an effective syringe pump designed to hold glass or plastic syringes of any make from 10 μl to 60ml
- Internal diameter of the syringe is used by the control program to calibrate the pump and deliver the volume and flowrate selected.
- Dispense Volume Mode: The pump keeps track of volume dispensed and automatically stops when target volume is reached.
- Run Mode: Pump runs at the set flow rate until stopped manually.
- Min. Flow Rate (10 µl syringe): 0.1 µl/hr
- Max. Flow Rate:

10 ml Syringe: 127 ml/hr 60 ml Syringe: 519 ml/hr

Tel: 916.429.7240 Fax: 916.429.2701